# cefar® peristim<sup>pro</sup>

|   |                             |     |
|---|-----------------------------|-----|
|  | USERS GUIDE .....           | 3   |
|  | MANUAL DEL USUARIO .....    | 25  |
|  | GEBRAUCHSANWEISUNG .....    | 47  |
|  | MODE D'EMPLOI .....         | 69  |
|  | ISTRUZIONI PER L'USO .....  | 91  |
|  | GEBRUIKERSHANDLEIDING ..... | 113 |
|  | BRUKSANVISNING .....        | 135 |

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## Urge incontinence

Urge incontinence means a sudden, strong urge to urinate followed by an immediate bladder contraction, resulting in an involuntary leakage of urine. Both men and women can be affected, particularly the elderly. One reason for this condition can be a disruption in the part of the nervous system that controls the bladder.

## Mixed incontinence

Mixed incontinence is a combination of stress incontinence and urge incontinence.

## Faecal incontinence

Faecal incontinence, also called anal or bowel incontinence, is the impaired ability to control passage of gas or stool. There are many possible causes of faecal incontinence, the most common is injury to the anal sphincter (ring-like muscle), for instance during childbirth or surgery, or damage to the nerves that control the anal sphincters. The condition usually becomes worse as people age.

## INCONTINENCE TREATMENT

Electrical stimulation via the pelvic nerves is a recognized treatment alternative for urinary incontinence. It is also proposed as a treatment method for faecal incontinence due to pelvic floor dysfunction or a poorly functioning anal sphincter.

When treating stress incontinence, the aim of the electrical stimulation is to mirror a voluntary muscle contraction and improve the function of the pelvic floor muscles. For urge incontinence, the aim is to inhibit involuntary bladder contractions by stimulating the nerves in the pelvic floor. When treating mixed incontinence, a stimulation appropriate for both urge incontinence and stress incontinence is used. For faecal incontinence, the aim is to improve bowel control by strengthening and toning pelvic floor muscles.

## 3. PRECAUTIONARY MEASURES

- Use only CEFAR-approved vaginal and anal probes intended for electrical stimulation.
- The stimulator may only be operated after instructions from medically trained personnel.
- The treatment is a personal ordination. Do not lend probes or stimulator to anyone else
- Only Cefar accessories should be used with the stimulator.

### WARNINGS

#### For internal use with a vaginal/anal probe

- People with implanted electronic equipment, such as pacemakers and intracardiac defibrillators, must not be treated with CEFAR PERISTIM PRO.
- Pregnant women must not be treated with the urology programs.
- People with extra-urethral incontinence (fistula, ectopic ureter) must not be treated with the urology programs.
- People with overflow incontinence due to outflow obstacle must not be treated with the urology programs.
- People with serious retention of urine in the upper urinary tract must not be treated with the urology programs.
- People with complete peripheral denervation of the pelvic floor must not be treated with the urology programs.
- Stimulation should not take place while the user is connected to high-frequency surgical equipment. It may cause burn injuries on the tissue, as well as problems with the stimulator.
- Do not use the stimulator in the vicinity of shortwave or microwave therapy equipment, since this may affect the output power of the stimulator.
- Keep the stimulator out of reach of children.

### CAUTION

#### For internal use with a vaginal/anal probe

- Patients with total/subtotal prolapsed uterus/vagina should be stimulated with greatest caution.
- Patients with urinary tract infections must be treated and clear of infection before starting therapy with the urology programs. Consult your doctor.
- If tissue irritation should occur, treatment should be temporarily discontinued. If problems continue, contact your health care provider. Hypersensitivity can occur in isolated cases. The problem usually disappears when probes or gel are changed to another type.
- Observe caution when using electrotherapy at the same time as the patient is connected to monitoring equipment with body worn electrodes. The stimulation might interfere with the signals to the monitoring equipment.
- Never open the battery cover during stimulation in order to avoid electrical shock.
- Turn off the stimulation or make sure that the amplitude for each channel is 0,0 mA before touching or removing the probes. Getting electrical stimulation through the fingers is unpleasant but not harmful.
- Observe caution when stimulating in the immediate vicinity of cellular phones that are switched on, since this may affect the output power of the stimulator.

#### For external use with surface electrodes

### 3.1. CONTRAINDICATIONS

- Implanted electronic devices. Do not use the device on patients who have a cardiac pacemaker, implanted defibrillator, or other implanted electronic device, because this may cause electric shock, burns, electrical interference, or death.
- TENS for Undiagnosed Pain. Do not use the device as a TENS device on patients whose pain syndromes are undiagnosed.

### 3.2. WARNINGS

- Consult with physician. Consult with the patient's physician before using the device, because the device may cause lethal rhythm disturbances to the heart in susceptible individuals.
- Skin condition. Apply stimulation only to normal, intact, clean, healthy skin.
- Long term effects. The long-term effects of chronic electrical stimulation are unknown.
- Stimulation location. Stimulation over Neck or Mouth. Do not apply stimulation over the patient's neck (especially the carotid sinus) or the patient's mouth, because this could cause severe muscle spasms resulting in closure of the airway, difficulty in breathing, or adverse effects on heart rhythm or blood pressure.
- Stimulation across Chest. Do not apply stimulation across the patient's chest, because the introduction of electrical current into the chest may cause rhythm disturbances to the patient's heart, which could be lethal.
- Across the head. Since the effects of stimulation of the brain are unknown, stimulation should not be applied across the head, and electrodes should not be placed on opposite sides of the head.
- Stimulation over compromised skin. Do not apply stimulation over open wounds or rashes, or over swollen, red, infected, or inflamed areas or skin eruptions (e.g., phlebitis, thrombophlebitis, varicose veins).
- Stimulation near cancerous lesions. Do not apply stimulation over, or in proximity to, cancerous lesions.
- Stimulation over eyes. Do not apply stimulation directly on the eyes.
- Environment. Electronic monitoring equipment. Do not apply stimulation in the presence of electronic monitoring equipment (e.g., cardiac monitors, ECG alarms), which may not operate properly when the electrical stimulation device is in use.
- Bath or Shower. Do not apply stimulation when the patient is in the bath or shower. Do not apply stimulation in humid atmosphere exceeding 75% of relative humidity. Sleeping. Do not apply stimulation while the patient is sleeping.
- Driving or operating machinery. Do not apply stimulation while the patient is driving, operating machinery, or during any activity in which electrical stimulation or involuntary muscle contraction can put the patient at risk of injury.
- Electrosurgical equipment or defibrillators. Disconnect the stimulation electrodes before using electrosurgical equipment or defibrillators. Otherwise skin burns may be caused below the electrodes and the device might be destroyed.
- Magnetic Resonance Imaging. Do not wear electrode or the device during Magnetic Resonance Imaging (MRI) scans as this may result in metal overheating and causing skin burns in the area of the electrode.
- Flammable or explosive environment. Do not use the device in areas where there is a risk of fire or explosion, such as oxygen-rich environments, in the vicinity of flammable anaesthetics, etc.
- Power supply. Never connect stimulation cables to an external power supply as there is a risk of electric shock.
- Near other equipment. Do not use the device beside or stacked on top of any other equipment. If you must use it side by side or on top of another system, you should check that the device works properly in the chosen configuration.
- Miscellaneous. Electrodes for Single Patient. Do not share electrodes with other persons. All users should have individual set of electrodes to prevent undesirable skin reactions or disease transmission.
- Accessories. Use this device only with the leads, electrodes, and accessories recommended by the manufacturer. Use of other accessories may adversely affect the performance of the device or may result in stronger electromagnetic emissions or reduce the electromagnetic immunity of the device.
- No Modification. No modification of the equipment is allowed.

### 3.3. PRECAUTIONS

- Supervision. Use this device only under the continued supervision of a licensed practitioner. Electrode placement and stimulation settings should be based on the guidance of the prescribing practitioner.


- **Manufacturer.** The manufacturer does not take any responsibility for any electrode placements other than recommended.
- **Pregnancy.** The safety of electrical stimulation during pregnancy has not been established.
- **Skin irritation.** Some patients may experience skin irritation or hypersensitivity due to the electrical stimulation or electrical conductive medium (gel). The irritation may be reduced by using an alternate conductive medium or alternate electrode placement. Some patients may experience redness under the electrodes after a session. This redness usually disappears within a few hours. Advise the patient to consult the clinician if the skin redness does not disappear after a few hours. Do not start another stimulation session in the same area if the redness is still visible. Don't scratch the redness area.
- **Heart Disease.** Patients with suspected or diagnosed heart disease should follow precautions recommended by their physicians.
- **Epilepsy.** Patients with suspected or diagnosed epilepsy should follow precautions recommended by their physicians.
- **Internal Bleeding.** Use caution when the patient has a tendency to bleed internally, such as following an injury or fracture.
- **After Surgery.** Use caution following recent surgical procedures when stimulation may disrupt the patient's healing process.
- **Over uterus.** Do not place electrodes directly over the uterus or connect pairs of electrodes across the abdomen if you are pregnant. The reason is that, theoretically, the current could affect the foetus's heart (although there are no reports of it being harmful).
- **Lack of sensation.** Use caution if stimulation is applied over areas of skin that lack normal sensation. Don't apply stimulation on patient unable to express themselves.
- **Stimulation.** The stimulator should only be used with skin electrodes intended for nerve and muscle stimulation. Muscle soreness may occur after stimulation but usually disappears within a week.
- **Hot casing or batteries.** Under extreme use conditions, some parts of the casing might reach up to 109 °F (43 °C). Use caution when manipulating the batteries right after device use or when holding the device. There is no particular health risk associated with this temperature besides your comfort.
- **Children.** Keep this device out of the reach of children.
- **Electrode Size.** Do not use electrodes with an active area less than 16 cm<sup>2</sup>, as there will be a risk of suffering a burn injury. Caution should always be exercised with current densities more than 2mA/ cm<sup>2</sup>.
- **Strangulation.** Do not wrap leadwires around your neck, and keep them out of the reach of children. Strangulation may result from entanglement in the leadwires.
- **Tripping.** Care should be used to avoid tripping on lead wires.
- **Damaged Device or Accessories.** Never use the device or any of its accessories if it is damaged (case, cables, etc.) or if the battery compartment is open as there is a risk of electric shock. Carefully inspect the lead wires and connectors prior to each use.
- **Inspect Electrodes.** Inspect electrodes before each use. Replace electrodes when they begin to deteriorate or lose adhesion. Poor contact between the electrodes and the patient's skin increases the risk of skin irritation or burns. Electrodes will last longer if used and stored according to instructions on electrode packaging. Attach the electrodes in such a way that their entire surface is in contact with the skin.
- **Foreign bodies.** Do not allow any foreign bodies (soil, water, metal, etc.) to penetrate the device and the battery compartment.
- **Batteries.** Do not carry batteries in a pocket, purse, or any other place where the terminals could become short-circuited (e.g. by way of paper clip). Intense heat could be generated and injury may result. Never open the battery cover during stimulation in order to avoid electrical shock. Remove the batteries from the device if you do not

# 4. OVERVIEW/CONTROL BUTTONS

## 1. ON/OFF

- Turns the stimulator on and off.
- Turns the stimulator off even when the key lock is activated. Can be used for terminating the stimulation at all times.

## 2. INCREASE (left and right channel)

- Increases the amplitude (intensity of stimulation). Press and hold the button to increase the amplitude continuously.

**Note!** Always increase the amplitude cautiously.

- Increases the number of minutes when setting the timer (right button).
- Used when turning the program lock on and off (left button)\*.
- Scrolls through choices in programming mode (left button)\*.

## 3. DECREASE (left and right channel)

- Decreases the amplitude (intensity of stimulation). Press and hold the button to decrease the amplitude continuously.

- Deactivates the key lock.
- Decreases the number of minutes when setting the timer (right button).
- Used when turning the program lock on and off (left button)\*.
- Scrolls through choices in programming mode (left button)\*.

## 4. PROGRAM

- Selects programs P1–P12.
- Pauses an ongoing program.
- Used when activating/deactivating the program lock\*.

## 5. TIMER

- Initiates the timer setting.

## 6. PROGRAMMING/CONFIRMATION\*

- Turns the stimulator into programming mode for custom programs P10–P12 when pressed for 2 seconds.
- Confirms settings in programming mode.

## 7. CABLE HOLDER

\*Professional use

# 7. PROGRAMS

## PRESET PROGRAMS

CEFAR PERISTIM PRO offers eight preset programs for stress, urge, mixed and faecal incontinence treatment and one preset program for relief of bladder-related pain.

Programs 1–7 should be used with a vaginal and/or anal probe and programs 8–9 with surface electrodes.

In programs 1–6 intermittent stimulation is used, which means that stimulation time alternates with rest time without stimulation: stimulation/rest/stimulation. The rest time allows the muscles to recover.

### Programs 1–3

#### STRESS INCONTINENCE

| Programnumber | Frequency | Pulse width | Stimulation time | Rest time | Recommended treatment time   |
|---------------|-----------|-------------|------------------|-----------|------------------------------|
| P1            | 50 Hz     | 200 µs      | 3 s              | 6 s       | 30 minutes, 3–5 times a week |
| P2            | 50 Hz     | 200 µs      | 5 s              | 10 s      | 30 minutes, 3–5 times a week |
| P3            | 50 Hz     | 200 µs      | 10 s             | 20 s      | 30 minutes, 3–5 times a week |

One cause of stress incontinence can be poorly functioning perineal muscles. Electrical stimulation of the perineal muscles with a vaginal and/or anal probe makes the muscles contract. The contractions help to identify and strengthen the weakened muscles. The stimulation must be as strong as possible without being painful, and it is helpful if the patient participates actively in the muscle contractions. The intensity should be sufficient enough to create a reflexive contraction of the anus. It is suggested that this training will be combined with the patient's own exercises for pelvic floor muscles.

### Programs 4–6

#### MIXED INCONTINENCE

| Programnumber | Frequency | Pulse width | Stimulation time | Rest time | Recommended treatment time   |
|---------------|-----------|-------------|------------------|-----------|------------------------------|
| P4            | 20 Hz     | 200 µs      | 3 s              | 6 s       | 30 minutes, 3–5 times a week |
| P5            | 20 Hz     | 200 µs      | 5 s              | 10 s      | 30 minutes, 3–5 times a week |
| P6            | 20 Hz     | 200 µs      | 10 s             | 20 s      | 30 minutes, 3–5 times a week |

Program 4–6 is used with a vaginal and/or anal probe for treating mixed incontinence. The stimulation should be as strong as possible without being painful. It is helpful if the patient participates actively in the muscle contractions.

## 9. CARE INSTRUCTIONS

Taking care of and cleaning the Cefar-Compex equipment is simple with the following instructions:

- For care and cleaning instructions of probes, see probe package.
- Self-adhesive multi-use electrodes are re-moistened if necessary with a few drops of water and kept air-tight (in a plastic bag) on protective paper when they are not in use.
- Rinse the carbon rubber electrodes and the skin with water after use. Do not use detergent for the electrodes.
- Never expose the stimulator to water. Wipe it off with a damp cloth if necessary.
- Do not jerk cables or connections.
- The cables are best preserved if left attached to the stimulator between sessions.
- No other maintenance of the device is required. Service life of the device can vary depending on usage conditions. Typical service life is 7 years.
- Never service the device while in use
- The device should be operated in temperatures between 10°C and 40°C, atmospheric pressures between 50 and 106 kPa, and relative humidity between 30% and 75%.
- The device should be transported and stored in temperatures between -40 °C and 70 °C, atmospheric pressures between 50 and 106 kPa and relative humidity between 10% and 90%.
- Contact manufacturer for assistance in setting up, using or maintaining the equipment or report events.

## 11. TECHNICAL DATA

CEFAR PERISTIM PRO is a stimulator with two non-independent channels intended for treatment of incontinence and bladder-related pain. It has nine preset programs and three custom programs.

Treatment with electrical stimulation requires the stimulation current to penetrate the resistance of the skin and the electrode, about 1000 ohms. CEFAR PERISTIM PRO can penetrate this resistance and maintain a current of up to 99.5 mA. With a change in load from 100 to 1000 ohms, the stimulation current changes less than 10% from the set value.

The stimulator operates on either two non-rechargeable 1.5 V AA batteries or two rechargeable 1.2 V AA batteries, recharged in a separate charger.

### CEFAR PERISTIM PRO

|   |   |
|---|---|
| Number of channels                                | 2 (non-independent)   |
| Constant current                                  | Up to a resistance of 1000 ohm<br>(increased load can reduce the maximum current) |
| Stimulation current/channel                       | 0–99,5 mA (charge maximum: 30µC)  |
| Waveform  | Symmetrical biphasic pulse, 100% compensated                                      |
| Number of preset programs                         | 9   |
| Number of custom programs                         | 3   |
| Stimulation forms                                 | Continuous stimulation<br>Intermittent stimulation                                |
| Max pulse duration                                | 300 µs  |
| Max frequency                                     | 120 Hz  |
| Timer   | 1 to 60 min/Off   |
| Environment for storage, use and shipping         | Temperature 10° C–40° C<br>Air humidity 30%–75%<br>Air pressure 700 hPa–1060 hPa  |
| Power source                                      | 2 x 1.5 V AA non-rechargeable or<br>2 x 1.2 V AA rechargeable batteries           |
| Current consumption for one channel, 80 Hz, 30 mA | 150 mA  |
| I r.m.s. max/channel                              | 27 mA   |
| Size  | 120 x 50 x 30 mm  |
| Weight  | ca. 180 g   |

## KEY TO THE SYMBOLS

Reference number, part number

Follow instruction for uses

Type BF applied parts

Council Directive 2012/19/EU concerning Waste Electrical and Electronic Equipment (WEEE) requires not to dispose of WEEE as municipal waste. Contact your local distributor for information of the unit and accessories

Keep the device dry

Minimum and maximum temperature indications to respect

This equipment complies with all requirements of the Medical Device Directive (93/42/EEC)

Manufacturer name and address, manufacturing year

Power/Pause

# 12. ELECTROMAGNETIC COMPATIBILITY (EMC) TABLES

| Guidance and Manufacturer's Declaration - Electromagnetic Emissions  |                                  |  |
|--|----------------------------------|--|
| The CEFAR PERISTIM PRO Electrotherapy System is intended for use in the electromagnetic environment defined below. The customer or the user of the CEFAR PERISTIM PRO Electrotherapy System should assure that it is used in such an environment |                                  |  |
| Emissions Tests  | Compliance                       | Electromagnetic Environment - Guidance   |
| RF emissions<br>CISPR 11   | Group 1                          | The CEFAR PERISTIM PRO Electrotherapy System uses RF energy only for its internal function. Therefore, its RF emissions are very low and are not likely to cause any interference in nearby electronic equipment.                                    |
| RF emissions<br>CISPR 11   | Class B                          | CEFAR PERISTIM PRO Electrotherapy System is suitable for use in all establishments, including domestic establishments and those directly connected to the public low voltage power supply network that supplies buildings used for domestic purposes |
| Harmonic emissions<br>IEC 61000-3-2  | Not Applicable - Battery powered |  |
| Voltage fluctuations<br>IEC 61000-3-3  | Not Applicable - Battery powered |  |

**Guidance and Manufacturer's Declaration - Electromagnetic Immunity**

The CEFAR PERISTIM PRO Electrotherapy System is intended for use in the electromagnetic environment specified below. The customer or the user of the CEFAR PERISTIM PRO Electrotherapy System should assure that it is used in such an environment

| Immunity Test  | IEC 60601 Test Level  | Compliance Level   | Electromagnetic Environment - Guidance  |
|--|---|--|---|
| Electrostatic discharge (ESD) IEC 61000-4-2  | ±6kV contact<br>±8kV air  | ±6kV contact<br>±8kV air   | Risk assessment on the CEFAR PERISTIM PRO Electrotherapy System indicates the compliance levels claimed are acceptable when ESD-precautionary measures are taken.   |
| Electrical fast transient/ burst IEC 61000-4-4   | ±2kV for power supply lines<br>±1kV for input/output lines  | Not Applicable - Battery powered<br><br>Not Applicable - signal lines less then 3 meters | Mains power quality should be that of a typical commercial or hospital environment.   |
| Surge IEC 61000-4-5  | + 1kV differential mode (line to line)<br>+ 2kV common mode (line to ground)  | Not Applicable - Battery powered   | Mains power quality should be that of a typical commercial or hospital environment.   |
| Voltage dips, short interruptions and voltage variations on power supply input lines<br><br>IEC 61000-4-11 | <5% $U_T$ (>95% dip in $U_T$ ) for 0,5 cycle<br>40% $U_T$ (60% dip in $U_T$ ) for 5 cycles<br>70% $U_T$ (30% dip in $U_T$ ) for 25 cycles<br><5% $U_T$ (>95% dip in $U_T$ ) for 5 sec | Not Applicable - Battery powered   | Mains power quality should be that of a typical commercial or hospital environment. If the user of the CEFAR PERISTIM PRO Electrotherapy System requires continued operation during power mains interruptions, it is recommended that the CEFAR PERISTIM PRO Electrotherapy System be powered from an uninterrupted power supply. |
| Power frequency (50/60Hz) magnetic field IEC 61000-4-8   | 3 A/m   | 3 A/m  | Power frequency magnetic fields should be at levels characteristic of a typical location in a typical commercial or hospital environment.   |

NOTE:  $U_T$  is the AC mains voltage prior to application of the test level.

### Guidance and Manufacturer's Declaration - Electromagnetic Immunity

The CEFAR PERISTIM PRO Electrotherapy System is intended for use in the electromagnetic environment defined below. The customer or the user of the CEFAR PERISTIM PRO Electrotherapy System should assure that it is used in such an environment.

| Immunity Test                 | IEC 60601 Test Level        | Compliance Level          | Electromagnetic Environment - Guidance   |
|-------------------------------|-----------------------------|---------------------------|--|
| Conducted RF<br>IEC 61000-4-6 | 3 Vrms<br>150 kHz to 80 MHz | [V1] V, where V1 = 3V     | <p>Portable and mobile RF communications equipment should be used no closer to any part of the CEFAR PERISTIM PRO Electrotherapy System, including cables, than the recommended separation distance calculated from the equation applicable to the frequency of the transmitter.</p> <p><b>Recommended separation distance</b></p> $d = \frac{3.5 \sqrt{P}}{[V1]}$   |
| Radiated RF<br>IEC 61000-4-3  | 3 V/m<br>80 MHz to 2,5 GHz  | [E1] V/m, where E1 = 3V/m | <p><math>d = \frac{3.5 \sqrt{P}}{[E1]}</math>    80 MHz to 800 MHz</p> <p><math>d = \frac{7 \sqrt{P}}{[E1]}</math>    800 MHz to 2,5 GHz</p> <p>where P is the maximum output power rating of the transmitter in watts (W) according to the transmitter manufacturer and d is the recommended separation distance in metres (m).</p> <p>Field strengths from fixed RF transmitters, as determined by an electromagnetic site survey,<sup>a</sup> should be less than the compliance level in each frequency range.<sup>b</sup></p> <p>Interference may occur in the vicinity of equipment marked with the following symbol:</p> <div style="text-align: center;">  </div> |

NOTE 1 At 80 MHz and 800 MHz, the higher frequency range applies.

NOTE 2 These guidelines may not apply in all situations. Electromagnetic propagation is affected by absorption and reflection from structures, objects and people.

<sup>a</sup> Field strengths from fixed transmitters, such as base stations for radio (cellular/cordless) telephones and land mobile radios, amateur radio, AM and FM radio broadcast and TV broadcast cannot be predicted theoretically with accuracy. To assess the electromagnetic environment due to fixed RF transmitters, an electromagnetic site survey should be considered. If the measured field strength in the location in which the CEFAR PERISTIM PRO Electrotherapy System is used exceeds the applicable RF compliance level above, the CEFAR PERISTIM PRO Electrotherapy System should be observed to verify normal operation. If abnormal performance is observed, additional measures may be necessary, such as reorienting or relocating the CEFAR PERISTIM PRO Electrotherapy System.

<sup>b</sup> Over the frequency range 150 kHz to 80 MHz, field strengths should be less than [V1] V/m.

**Recommended separation distances between portable and mobile RF communications equipment and the CEFAR PERISTIM PRO Electrotherapy System**

The CEFAR PERISTIM PRO Electrotherapy System is intended for use in an electromagnetic environment in which radiated RF disturbances are controlled. The customer or the user of the CEFAR PERISTIM PRO Electrotherapy System can help prevent electromagnetic interference by maintaining a minimum distance between portable and mobile RF communications equipment (transmitters) and the CEFAR PERISTIM PRO Electrotherapy System as recommended below, according to the maximum output power of the communications equipment.

| Rated maximum output power of transmitter P (W) | Separation distance according to frequency of transmitter d (m)       |   |   |
|---|---|---|---|
|   | 150 kHz to 80 MHz<br>$d = \frac{35}{[V]} \sqrt{P}$<br>(where V1 = 3V) | 80 MHz to 800 MHz<br>$d = \frac{35}{[E]} \sqrt{P}$<br>(where E1 = 3V/m) | 800 MHz to 2,5 GHz<br>$d = \frac{7}{[E]} \sqrt{P}$<br>(where E1 = 3V/m) |
| 0,01  | 0,12  | 0,12  | 0,23  |
| 0,1   | 0,38  | 0,38  | 0,73  |
| 1   | 1,2   | 1,2   | 2,3   |
| 10  | 3,8   | 3,8   | 7,3   |
| 100   | 12  | 12  | 23  |

For transmitters rated at a maximum output power not listed above, the recommended separation distance d in metres (m) can be estimated using the equation applicable to the frequency of the transmitter, where P is the maximum output power rating of the transmitter in watts (W) according to the transmitter manufacturer.

NOTE 1 At 80 MHz and 800 MHz, the separation distance for the higher frequency range applies

NOTE 2 These guidelines may not apply in all situations. Electromagnetic propagation is affected by absorption and reflection from structures, objects and people.

**FCC REQUIREMENTS**

**Part 15 of the FCC Requirements**

|   |   |
|---|---|
| This device complies with Part 15 of the FCC Rules. Operation is subject to the following 2 conditions: | <ul style="list-style-type: none"> <li>This device may not cause harmful interference</li> <li>This device must accept any interference received, including the interference that may cause undesired operation.</li> </ul> |
| FCC ID  | T9J-RN42  |
| Contains Transmitter Module IC  | 6514A-RN42  |

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## 13. LITERATURVERZEICHNIS

### **Fall M, Lindström S**

Electrical stimulation: A physiologic approach to the treatment of urinary incontinence. *Journal of Urologic Clinics of North America*; 2, Vol 18: 393-407, 1991

**Fall et al.** Transcutaneous Electrical Nerve Stimulation in classic an nonulcer interstitian cystitis. *Urologic clinics of North America Vol 21, No 1, Feb 1994*

**Walsh et al.** Non-invasive antidromic neurostimulation, a simple effective method for improving bladder storage. *NeuroUrol Urodyn 2001;20(1):73-84*

**Amarenco et al.** Urodynamic effect of acute transcutaneous posterior tibial nerve stimulation in overactive bladder. *J Urol.2003 Jun;169(6):2210-5*

### **Zöllner-Nielsen M, Samuelsson S.M**

Maximal electrical stimulation of patients with frequency, urgency and urge incontinence. Report of 38 cases. *Acta Obstet Gynecol Scand 1992;71:629-631*

### **Sand P.K, Richardson D.A Et al.**

Pelvic floor electrical stimulation in the treatment of genuine stress incontinence: A multicenter, placebo-controlled trial. *Am J Obstet Gynecol Vol 173, No 1, 1995*

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